

OMNigene™•SALIVA DNA and RNA (OMR-610) purification protocol for viral total nucleic acids using Applied Biosystems MagMAX™ Viral/Pathogen Nucleic Acid Isolation Kit

This laboratory protocol is used for the preparation of a saliva sample collected and stabilized in OMNigene™•SALIVA DNA and RNA (OMR-610) for subsequent extraction of viral total nucleic acids using MagMAX™ Viral/Pathogen Nucleic Acid Isolation Kit (Applied Biosystems, Cat. No. A42352).

Equipment and reagents

- MagMAX Viral/Pathogen Nucleic Acid Isolation Kit (Applied Biosystems)
- Equipment and reagents listed are referenced in the MagMAX Viral/Pathogen Nucleic Acid Isolation Kit (Applied Biosystems) user guide.
For a detailed extraction protocol, refer to the user guide included in the extraction kit.

Procedure

Sample preparation steps
1. Shake OMR-610 samples vigorously for 10 seconds once received in the laboratory.
2. Incubate the entire sample in the original vial at 50°C for 1 hour in a water bath or for 2 hours in an air incubator. Note: Incubation may be performed any time between sample receipt and extraction. This step does not need to be repeated for extraction of subsequent aliquots.
3. Continue with the MagMAX Viral/Pathogen Nucleic Acid Isolation Kit (Applied Biosystems) instructions, starting at step 1.

Technical support is available Monday to Friday (9:00 a.m.-5:00 p.m. ET):

- Toll-free (North America): 1.866.813.6354, option 6
- All other countries: 613.723.5757, option 6
- Email: support@dnagenotek.com

Some DNA Genotek products may not be available in all geographic regions.
OMNigene-SALIVA DNA and RNA (OMR-610) is for Research Use Only, not for use in diagnostic procedures.
OMNigene and DNA Genotek are trademarks of DNA Genotek Inc. All other brands and names contained herein are the property of their respective owners.
All DNA Genotek protocols, white papers and application notes are available in the support section of our website at www.dnagenotek.com.

Not available for clinical diagnostic use in the United States.