

For Animal Use Only

Question	Possible cause	Explanation/solution
Fluorescence and absorbance quantification values do not match.	RNA is present in sample.	It is normal to obtain different values when measuring yields with absorbance vs fluorescence, as absorbance will measure both DNA and RNA content, while fluorescence will only measure DNA content. We suggest quantifying by fluorescence as it is a more accurate way of measuring DNA concentrations. If quantification by absorbance is necessary, we recommend RNase-treating the samples. A detailed protocol is described in PD-PR-040, <i>RNA removal by double-RNase digestion</i> , and can be found on our website.
Low OD ratios.	The PG-L2P purifier was not added.	The PG-L2P purifier must be added to precipitate impurities and inhibitors from the sample.
	Incorrect blank.	Ensure that you blank your spectrophotometer using the same buffer that was used to dissolve the DNA pellet (e.g., TE).
	Absorbance readings were not corrected for turbidity.	When measuring by absorbance it is essential to subtract the A_{320} reading from your A_{280} and A_{260} readings before calculating the 260/280 ratio. This eliminates the contribution of any turbid material present in your sample. Note that this turbid material is inert and should not affect your downstream applications. The formula for this subtraction can be found in the purification protocol. Some spectrophotometers may perform this correction automatically. Please see the documentation for your instrument.
	Residual ethanol not removed.	Residual ethanol should be collected at the bottom of the tube by brief centrifugation and removed using a pipette. If residual ethanol is evaporated, impurities left in the wash will concentrate and be left behind in the sample affecting the 260/280 ratio.
	Post PG-L2P purifier pellet was disturbed.	The pellet contains turbid impurities. Even a small carryover of this pellet can reduce ratios. If accidentally disturbed, the tube should be re-centrifuged before removing the supernatant. A longer period of centrifugation, (up to 15 minutes) after the addition of PG-L2P and 10 minutes incubation in ice may be beneficial in reducing the turbidity (high A_{320}) of the final DNA solution.
	Turbid material in extracted DNA.	There are several methods to help remove insoluble particulate from the DNA. These include high-speed centrifugation and chloroform treatment. More details are available in MK-AN-017, <i>From turbidity to clarity: Simple methods to improve the A_{260}/A_{280} ratio of Oragene-purified DNA samples</i> , which can be found on our website.
Storage of raw PG-100 samples and extracted DNA.	Short-term storage of raw sample.	Unpurified samples may be stored at room temperature in their original tube for up to 1 year.
	Long-term storage of raw sample.	For long-term storage, we recommend transferring the samples to a cryo-tube with an O-ring and storing at -20°C or -80°C . Samples may be split into aliquots and stored in microcentrifuge tubes. To ensure sample homogeneity, heat the entire sample for 2 hours at 50°C before aliquoting.
	Storage of purified DNA sample.	Store purified DNA in 1×TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 8.0 or greater). Samples may be stored short-term (weeks) at 4°C or long term frozen at -20°C .

