


## Storage and aliquoting recommendations for samples collected in OMNIGENE™·GUT Dx (OMD-200)

This document contains the recommendations on how to store OMNIGENE™·GUT Dx collection kits after sample collection and how to aliquot samples for storage or processing.

Additional information regarding this product can be found in the OMNIGENE™·GUT Dx (OMD-200) product handbook, PD-HB-00023.

### OMNIGENE™·GUT Dx (OMD-200) post-collection storage recommendations

	Room temperature storage (15°C / 25°C/59°F / 77°F)	Freezer storage (-20°C or -80°C/-4°F or -112°F)
<b>OMNIGENE™·GUT Dx</b> (OMD-200) 	<b>Storage at room temperature:</b> Fecal samples collected with OMNIGENE™·GUT Dx (OMD-200) devices can be stored at room temperature for up to <b>30 days</b> .  <b>IMPORTANT:</b> Storing at 4°C (39°F) is NOT recommended for fecal samples collected with OMNIGENE™·GUT Dx devices.	<b>Storage at -20°C or -80°C (-4°F or -112°F):</b> Nucleic acids should be extracted within 30 days of the sample being collected. Extracted samples should be stored at -20°C or -80°C (-4°F or -112°F) as per your standard operating procedures.
		<b>Freeze-thaw cycles:</b> OMNIGENE™·GUT Dx will preserve the microbial community structure and maintain DNA integrity for up to <b>5 freeze-thaw cycles</b> . <sup>1</sup>

### OMNIGENE™·GUT Dx (OMD-200) aliquoting recommendations





Recommendations related to aliquoting fecal (stool) samples collected and stabilized in OMNIGENE™·GUT Dx devices.

#### Equipment required


- Pipettors and wide-bore pipettor tips (e.g., VWR 89049-160)
- Bench top vortex
- Sterile cryovials with O-rings
  - The final collected volume of an OMD-200 tube is greater than 1.8 mL (maximum volume that should be placed in a single 2 mL cryovial)
  - Popular practice is to divide homogenized fecal sample evenly between 2 x 2 mL cryovials with O-rings.

*Continued on next page*

## Procedure

	<ol style="list-style-type: none"> <li>1. Collect a fecal sample according to the OMNigene™•GUT Dx instructions. (DNA Genotek, OMD-200 Instructions for Use, PD-PR-01261.)</li> </ol>
	<ol style="list-style-type: none"> <li>2. Vortex the sample vigorously for 60 seconds (medium setting). This action typically will break up the matrix, making it more liquid and visibly homogenous.</li> </ol>
	<ol style="list-style-type: none"> <li>3. With the purple cap still screwed on, unscrew the yellow portion of the tube top and set aside on a clean surface.</li> </ol>
	<ol style="list-style-type: none"> <li>4. Using a wide-bore 1000 µL pipette tip, pipette out the sample and transfer it to a cryovial with an O-ring. If there is no wide-bore pipette tip, cut off the end of a pipette tip after sterilizing the scissors.</li> </ol>

## Troubleshooting

	<p>If the sample is still too viscous to pipette, put the purple cap and yellow portion back on the tube tightly. Contact DNA Genotek to obtain OMNigene™ Liquefaction Reagent (OM-LQR 400/1600) and follow instructions 1 through 5 as outlined in the OMNigene™•GUT Dx collection device bacterial DNA purification protocol (PD-PR-00968, DNA Genotek).</p>
---	--

For the optimized extraction protocol, see: OMNigene™•GUT Dx Collection device bacterial DNA purification protocol using QIAGEN QIAamp PowerFecal Pro DNA kit (PD-PR-00968, DNA Genotek).

1 DNA Genotek, OMNigene•GUT Dx (OMD-200) product handbook, PD-HB-00023.

**Technical support is available Monday to Friday (9h00 to 17h00 ET):**

Toll-free (North America): 1.866.813.6354, option 6

All other countries: +1.613.723.5757, option 6

Email: [support@dnagenotek.com](mailto:support@dnagenotek.com)