



		Pre-collection			Post-collection		
Catalogue numbers	Saliva volume	Package dimensions and contents — WEIGHT 34-38g [†]			Total sample volume ^{††}	Tube dimensions (mm)	
		Clamshell (mm)	Components	Collector			
Without 2D bottom barcode	OG(X)-600		<p>Instructions for use</p> <p>Cap</p>		4 mL		
	OG(X)-675				.75 mL via sponge		1.5 mL
	OG(X)-610				1 mL		2 mL
	ON-600				.5 mL		
With 2D bottom barcode	OG(X)-600.S2P		<p>Sponges for assisted collection (OG(X)-675 products only)</p>	<p>OR</p> <p>2D bottom barcode</p>	4 mL		
	OG(X)-675.S2P				.75 mL via sponge		1.5 mL
	OG(X)-610.S2P				1 mL		2 mL
	ON-600.S2P				.5 mL		

† Weight is inclusive of all product configurations within referenced catalogue numbers.

†† Total sample volume includes saliva and Oragene stabilization solution. May vary based on volume of saliva collected.

© 2017 DNA Genotek Inc., a subsidiary of OraSure Technologies, Inc., all rights reserved. Patent (www.dnagenotek.com/legalnotices)