


Oral sample collection
EASY • PAINLESS • PROVEN

Collect superior samples for your genetic analysis

ORAc collect[®] RNA is an all-in-one system for the collection, stabilization and transportation of RNA from human saliva for *In vitro* diagnostic use. The user-friendly kit design makes it ideal for self or assisted collection from donors of any age, in a clinic or at-home. ORAc collect RNA stabilizes RNA at room temperature for months and is compatible with downstream expression analysis applications such as RT-qPCR, microarray, and RNAseq.

Applicable for:

- miRNA Diagnostics
- Human mRNA and miRNA expression profiling
- Developing “expression signatures” of disease
- Researching non-human forms of RNA found in saliva (i.e., viral and bacterial RNA)

Specifications:

- Catalogue number: ORE-100
- Sample stability: 60 days
- Shelf-life: 12 months
- Tube diameter: 15 mm
- Tube height (uncapped): 92 mm
- Capped tube height (post collection): 103 mm
- Pre-collection kit weight: 7.05 g


Superior samples • Proven performance

Some DNA Genotek[™] products may not be available in all geographic regions.

*ORAc collect is a registered trademark and DNA Genotek is a trademark of DNA Genotek Inc. All other brands and names contained herein are the property of their respective owners. Patent (www.dnagenotek.com/legalnotices) • © 2020 DNA Genotek Inc., a subsidiary of OraSure Technologies, Inc., all rights reserved. PD-BR-00365 Issue 2/2020-09


For In Vitro Diagnostic Use

DNA Genotek Inc.
3000 - 500 Palladium Drive
Ottawa, ON, Canada K2V 1C2
Subsidiary of OraSure Technologies, Inc.

Toll-free (North America): 1.866.813.6354
Tel.: +1.613.723.5757 • Fax: +1.613.723.5057
info@dnagenotek.com
www.dnagenotek.com

DNAGENOTEK[™]