

Product change notification

Product affected:  **ORACollect•DNA**

Catalog number: OCR-100

Description: ORACollect•DNA

Effective date of change: August 2017

Change description:

ORACollect•DNA storage instructions are changing from 15°C – 30°C to 15°C – 25°C. Additionally, the product shelf-life is changing from 30 months to 24 months from date of manufacture.

Reason for change:

DNA Genotek is aligning storage instructions and product shelf-life across all ORACollect devices.

Customer action:

No action required. ORACollect•DNA products manufactured before August 2017 are not affected and can be used per the device labelling.

Please contact info@dnagenotek.com if you have any issues or concerns pertaining to this product change notification.