



Sample storage box CONVENIENT • DURABLE • REUSABLE

For easy and convenient sample storage

DNA Genotek's sample storage boxes are designed to provide convenient sample storage. Supplied flat-packed prior to use, the storage boxes save valuable laboratory shelf-space and eliminate the hassle of stacking. Made of a durable moisture resistant material, they outlast cardboard alternatives making them ideal for re-use. Each storage box holds 36 Oragene®, ORAcollect®, OMNIgene® or Performagene™ collection tubes.

- **Flat-pack reduces storage space**
- **Fits in standard freezer racks**
- **Moisture resistant material ideal for re-use**

Product specifications

Suitable for: 15 mL tubes
 Storage capacity: 36 tubes (19.8 mm diameter per well)
 Dimensions: 14.5 cm length × 14.7 cm width × 12 cm height
 Packaging: 10 storage boxes per unit, supplied complete with dividers



Catalog #: SO-1

For more information contact
info@dnagenotek.com

Some DNA Genotek products may not be available in all geographic regions.
 *Oragene, ORAcollect and OMNIgene are registered trademarks and Performagene is a trademark of DNA Genotek Inc.
 All other brands and names contained herein are the property of their respective owners.

*Superior samples
 Proven performance*