



Product change notification

Products affected:  **PERFORMAGENE™**

Catalog number: PG-100

Description: Performagene for collection of animal DNA samples

Effective date of change: May 1, 2019

Customer action: No required action is expected. Should you use the YYMMLL or SSSS (serial number) information embedded within the barcode, please take note of the shift in location within the new barcode nomenclature.

Reason for change: DNA Genotek is implementing manufacturing efficiencies to allow an increase in manufactured lot size.

Change description: The 14 character barcode is a unique identifier which contains data related to the manufacturing of the device, specifically product code, product SKU, device lot and serial number.


The 4 digit serial number will be changing to a 5 digit serial number, and the 3 digit product code will be reduced to 2 digits. There is no change to the overall number of digits (14) in the barcode.

Current 14 character nomenclature



3 Product code (ie., Performagene)
100 **Product SKU (PG-100)**
YYMMLL Lot # (year, month, lot#)
SSSS **4 digit unique serial number**

New 14 character nomenclature



3 Product code (ie., Performagene)
10 **Product SKU (ie., 10 for PG-100)**
YYMMLL Lot # (year, month, lot#)
SSSSS **5 digit unique serial number**

Note: pink denotes variable content.

Please contact info@dnagenotek.com if you have any issues or concerns pertaining to this product change notification.