


**Oral sample collection**
**EASY • PAINLESS • PROVEN**

## Collect superior samples for your genetic analysis

All expression analysis starts with the collection of RNA samples. ORAc collect<sup>®</sup>RNA (OR-100) is a Research Use Only, all-in-one system for the collection, stabilization and transportation of RNA from saliva.


- Quick and reliable collection of RNA
- Easy collection makes it ideal for self or assisted collection
- Stabilizes human, viral and bacterial RNA found in saliva at room temperature
- Compatible with downstream applications (e.g., RT-qPCR, microarray etc.)

### Applicable for:

- Human mRNA expression profiling
- Developing “expression signatures” of disease
- Validating detection of infections (e.g., influenza, coronavirus, etc.)

### Specifications:

- Catalogue number: OR-100
- Sample stability: 60 days
- Shelf-life: 12 months


*Superior samples • Proven performance*

Some DNA Genotek<sup>™</sup> products may not be available in all geographic regions.

<sup>®</sup>ORAc collect is a registered trademark of DNA Genotek Inc. All other brands and names contained herein are the property of their respective owners.

Patent ([www.dnagenotek.com/legalnotices](http://www.dnagenotek.com/legalnotices)) • © 2020 DNA Genotek Inc., a subsidiary of OraSure Technologies, Inc., all rights reserved. PD-BR-00364 Issue 1/2020-03

**For Research Use Only**  
**Not for use in diagnostic procedures**  
**For collection of human RNA**

DNA Genotek Inc.  
 3000 - 500 Palladium Drive  
 Ottawa, ON, Canada K2V 1C2  
 Subsidiary of OraSure Technologies, Inc.

Toll-free (North America): 1.866.813.6354  
 Tel.: +1.613.723.5757 • Fax: +1.613.723.5057  
[info@dnagenotek.com](mailto:info@dnagenotek.com)  
[www.dnagenotek.com](http://www.dnagenotek.com)

**DNAGENOTEK<sup>™</sup>**