

Product change notification

Products affected:  **ORAgene™•Dx**

Catalog numbers/SKUs: OGD-500, OGD-510, OGD-575, OGD-600, OGD-610, OGD-675

Description: Oragene™•Dx

Effective date of change: Effective immediately

Change description:

As part of the continuous improvement process at DNA Genotek, the post-collection DNA stability claim for the Oragene™•Dx family of devices has been extended from 1 year to 5 years, providing additional sample storage flexibility at room temperature for customers aiming to reduce their carbon footprint.

This change is applicable globally, with the exception of the United States, and is retroactive for Oragene™•Dx devices in the field. Previously collected saliva samples can now be stored for up to 5 years at room temperature.

Customer action:

No action is required.

For any questions pertaining to this product change notification, contact info@dnagenotek.com.

