

Product change notification

Products affected: Bio-specimen bag


Catalog numbers/SKUs: MO-5

Change description: The MO-5 bio-specimen bag is a specialized leak-proof bag containing absorbent material. The bag is intended to support the transportation of collected bio-specimen samples. The MO-5 bio-specimen bag features an adhesive zone used to secure the bio-specimen sample inside. The adhesive zone is being reduced from a 30-mm zone to a 20-mm zone.

This change does not affect the bio-specimen bag's leak-proof performance or its intended use.

It is possible that you may receive MO-5 bio-specimen bags with the 30-mm adhesive zone until the inventory is depleted.

Current adhesive zone



New adhesive zone



Customer action: No action is required.

For any issues or concerns pertaining to this product change notification, email info@dnagenotek.com.