

DNA from saliva

Collect high-quality DNA samples for molecular diagnostics with an FDA-cleared device.

Genetic testing starts with the collection of DNA samples. Use the Oragene™•Dx device for the proven collection, stabilization and transportation of DNA from saliva*. For donors who cannot comply with blood collection or who are unable to spit, the Oragene™•Dx (OGD-675) device allows for the minimally invasive collection of DNA from saliva using absorbent sponges.




- **Painless and reliable DNA collection in any setting – at home and in a clinic**
- **High-quality, high-quantity DNA compatible for a range of downstream applications[†]**
- **No cold chain required – transport DNA samples at ambient temperature via standard postal system**
- **Integrated 1D barcode for laboratory workflow efficiency**



Oragene™•Dx devices[†] have an FDA 510(k) general clearance for prescription use (Rx) and direct-to-consumer (DTC) uses. Molecular diagnostic manufacturers could reduce the time and regulatory burden of their own FDA submissions by leveraging the 510(k) clearances of Oragene™•Dx devices.

Robustly validated for technical performance and successful use by naive participants



Supplemental services

- Comprehensive custom offerings, including branding and fulfillment services
- Regulatory guidance and technical support
- Access instructions for use in multiple languages

Performance data

Attributes	SKU		
	OGD-600	OGD-610	OGD-675
Collection method	Self-collection	Self-collection	Assisted collection
Saliva volume (mL)	2	1	0.75
Mean DNA yield (μg)**	58.5	36.0	13.5
Molecular weight	> 23 kb	> 23 kb	> 23 kb
Mean concentration (ng/μL)	68.1	76.3	41.1
A ₂₆₀ /280 ratio	1.7	1.7	1.7

Device contents and packaging specifications


Shelf life: 30 months (15°C to 30°C / 59°F to 86°F)***

Full product specifications are available at www.dnagenotek.com. [PD-BR-233].








Mailing accessories available.



Self-collection (OGD-600/OGD-610)



Assisted collection (OGD-675)



For full collection instructions, visit www.dnagenotek.com

*For Prescription Use [k152464] and Over-the-Counter (Direct-to-Consumer) Use [k212745].

**DNA Genotek. OrAgene™-Dx product handbook. [PD-HB-1].

*** Shelf life is 30 months from date of manufacture at specified storage conditions.

Some DNA Genotek products may not be available in all geographic regions.

OrAgene-Dx is a trademark of DNA Genotek Inc., registered in various jurisdictions. All other brands and names contained herein are the property of their respective owners. Patent (www.dnagenotek.com/legalnotices)